



For Professional Scar Therapy

Silicone Sheet

Silicone Gel

- Affects both old and new scars
- Softens and flattens raised scars
- Relieves the itching and discomfort associated with scars
- Reduces redness and discoloration
- Prevents abnormal and excessive scar formation



Under License of Eska Pharma GmbH Germany



Class I Medical Device



Scar Gel

Medical Silicone Gel

DG DermTM Scar Gel

DG DermTM Silicone Scar Gel is a topical liquid silicone gel, made with a proprietary formula designed to dry quickly and leave a thin, transparent film that does not flake off on the skin. Using Scar Gel is an easy, safe, and convenient way to improve abnormal scars without the need to go through painful invasive procedures.

Features & Benefits

- Prevents abnormal scarring
- Quick Absorption
- Odorless, Colorless
- Reduces itching and discomfort associated with scars
- Suitable for all skin types and age
- Clinically tested to work on both old and new scars
- Easy application

How to Use

Apply a thin layer of Scar Gel on clean and dry scar. Use it 2-3 times daily for 2 months. Older scars may take longer.

Note: Do not use on open wounds. Keep out of reach of children.



• Keloid scars



• Hypertrophic Scars



• Burn Scars



• Normotrophic Scars



• Surgical Scars



• Traumatic Scars



• Acne Scars



• Stretch Marks

Ingredients

Medical grade silicone gel, Carbomer, Emulsifier

Side Effects

Due to the fact that silicone products do not have systemic absorption, they rarely cause secondary side effects.



Advanced Silicone Scar Sheet

Clinically proven to help fade, flatten & soften scars

DG Derm[®] Scar Sheet

DG Derm[®] Silicone Scar Sheet is a transparent, soft, and elastic sheet, mainly made of medical grade silicone. The sheet is effective in preventing the formation of new scars and reducing the appearance of existing ones. Silicone hydrates the scar tissue and regulates the production of collagen fibers. It has also been shown to reduce redness and make the scar flatter, softer, and paler.

DG Derm Silicone Scar Sheet is suitable for keloids, hypertrophics, burns, surgical scars and traumatic scars.

Features & Benefits

- 100% Medical grade silicone
- Various in shape and size
- Durable and elastic
- Comfortable and reusable
- Optimal scar healing by providing elastic protection
- No irritation
- Latex free
- Regain beauty

Mechanism of Action




How to Use

The silicone sheet is reusable and should be applied on clean and dry skin. Remove the plastic backing from the silicone sheet and put the sticky side over the scar. The sheet can be fixed in place by hypoallergenic cloth tape in case of losing stickiness over time.


The sheet can be washed with mild soap and water after taking it off if dirty.

The minimum effective usage duration is 12 hours daily, and expendable up to 20 hours. The treatment protocol is 8 weeks on the condition of constant application.

Note: Do not use on open wounds.




• C-section



• 2 months later



• Keloids



• 2 months later

The scar becomes softer and close to the surrounding skin color.



Silicone sheet shape	Silicone sheet size (cm)
	5 x 2/5 cm
	7 x 3 cm
	12 x 3/5 cm
	15 x 5 cm
	30 x 5 cm
	12 x 10 cm
	21 x 12 cm



No. 3, Dena Building, Farhang Shahr St., Shiraz, Iran
 Tel : +98 71-91002930 | Fax: +98 71-36302771
 www.dgdena.com | info@dgdena.com

To Know More

