

Dena Scrub™ Surgical Scrub Brush



Features :

- Available in dry or impregnated with cleansing solutions
- Supplied in handy upright dispenser box
- Polythene peel pouch, easy-open for dry or wet hands
- Flexible soft nylon bristles
- Polyethylene nail cleaner
- Plush polyester foam sponge



10993-10:2010



13485:2016

Product Description :

D.G.Dena offers a portfolio of high-quality Surgical Scrub Brushes for optimal and secure surgical hand preparation.

Our Scrub Brush helps reduce the risk of irritation or skin damage, thanks to its soft bristles. D.G.Dena offers dry sterile and impregnated Scrub Brushes, which evenly dispense the antiseptic solution.



Product Model	Product Description	Quantity per Box
D001	Dry, Solution-free and Sterile	30 pieces
D002	Impregnated with Chlorhexidine Digluconate 4%	30 pieces
D003	Impregnated with Povidone Iodine 7.5%	30 pieces

Our Portfolio of Scrub Brushes Offers :

- Easy-open packaging for wet or dry hands
- Dry sterile or impregnated with Chlorhexidine Digluconate (CHG) and Povidone Iodine (PVP-I)
- Easy-to-use wall-mounted dispensers.
- Incorporation of popular and efficacious antimicrobial formulations, including CHG and PVP-I
- Proprietary Polyethylene soft bristles to enhance mechanical action while reducing skin irritation
- Soft sponge to provide gentle cleaning on irritated areas or when firm mechanical action is not required
- Ergonomic brush design for greater user comfort and effectiveness

